MT 200

IB MT 200 S-V11 2420, Revision Date: 2020-06-03

Description of this Thermometer

- ON/OFF button
- Display
- Battery compartment cover
- Measuring sensor

Important Safety Instructions

- Follow instructions for use. This document provides important product operation and safety information regarding this device. Please read this document thoroughly before using the device and keep for future reference. This device may be used only for measuring body temperature!
- The minimum measurement time until the beep is heard must be maintained without exception!
- Ensure that children do not use this device unsupervised; some parts are small enough to be swallowed.
- Do not use this device close to strong electromagnetic fields such as mobile telephones or radio installations. Keep a minimum distance of 3.3 m from such devices when using this device.
- Protect the device from impact and dropping!
- Avoid bending the tip more than 45°!
- Avoid ambient temperatures above 60 °C. NEVER boil this device!
- Use only the commercial disinfectants listed in the section «Cleaning and Disinfecting» to clean the instrument. The instrument must be intact when immersed in liquid disinfectant.
- We recommend this device is tested for accuracy every two years or after mechanical impact (e.g. being dropped). Please contact Micro-life-Service to arrange the test.

⚠ WARNING: The measurement result given by this device is not a diagnosis! Do not rely on the measurement result only.

Batteries and electronic devices must be disposed of in accordance with the locally applicable regulations, not with domestic waste.

Read the instructions carefully before using this device.

Type BF applied part

Turning on the Thermometer

To turn on the thermometer, press the ON/OFF button ①; a short beep signals «thermometer ON». A display test is performed. All segments should be displayed.

The last measurement reading will be shown on the display ② automatically for 2 seconds with the «M» icon.

Then at an ambient temperature of less than 32 °C, an «L» and a flashing «°C» appear at the display field ②. The thermometer is now ready for use.

Function Test

Correct functioning of the thermometer is tested automatically each time it is turned on. If a malfunction is detected (measurement inaccuracy), this is indicated by «ERR» on the display, and a measurement becomes impossible. In this case, the thermometer must be replaced.

Using the Thermometer

Choose the preferred measuring method. When taking a measurement, the current temperature is continuously displayed and the «°C» symbol flashes. If the beep is heard 10 times and the «°C» is no longer flashing, the predictive end-temperature has been determined and the thermometer can be read now.

10 short beeps will sound when the temperature is higher than 37.5 °C in order to alert the patient that he/she may have fever. Reference: Oral temperature.

To achieve comparable results allow a 1 minute interval time between measurements.

To prolong the battery life, turn off the thermometer by briefly pressing the ON/OFF button ①. Otherwise the thermometer will automatically turn off after about 10 minutes.

Measuring methods / Normal body temperature

► In the mouth (oral) / 35.5 - 37.5 °C

Position the thermometer in one of the two pockets under the tongue, to the left or right of the root of the tongue. The measuring sensor ④ must be in good contact with the tissue. Close your mouth and breathe evenly through the nose to prevent the measurement from being influenced by inhaled/exhaled air.

Approx. measuring time: 10 seconds!

► **In the anus (rectal) / 36.6 - 38.0 °C**
Carefully insert the measuring sensor ④ of the thermometer 2 to 3 cm into the anal aperture.

Approx. measuring time: 10 seconds!

► **In the armpit (axillary) / 34.7 - 37.3 °C**

To receive more reliable results we recommend measuring temperature orally or rectally. A minimum measurement time of 5-10 minutes is recommended regardless of the beep sound.

Cleaning and Disinfecting

Name: Isopropyl alcohol 70 %; immerse: max. 24 hours.

Battery Replacement

When the «▼» symbol (upside-down triangle) appears in the display, the battery is flat and needs replacing. To replace the battery remove the battery compartment cover ③ from the thermometer. Insert the new battery with the +» at the top. Make sure you have a battery of the same type to hand. Batteries can be purchased at any electrical store.

Technical Specifications	
Type:	Predictive maximum thermometer
Measurement range:	32.0 °C to 42.9 °C <p>Temp. < 32.0 °C: display «L» for low (too low) Temp. > 42.9 °C: display «H» for high (too high)</p>
Measurement accuracy:	± 0.1 °C between 34 °C and 42 °C
Operating conditions:	10 - 40 °C; 15-95 % relative maximum humidity
Storage conditions:	-25 - +60 °C; 15-95 % relative maximum humidity
Battery:	1.5/1.55 V; SR41
Battery lifetime:	approx. 2700 measurements (using a new battery)
IP Class:	IP67
Reference to standards:	EN 12470-3, clinical thermometers; ASTM E1112; IEC 60601-1; IEC 60601-1-2 (EMC); IEC 60601-1-11
Expected service life:	5 years or 10000 measurements

This device complies with the requirements of the Medical Device Directive 93/42/EEC.

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